




Faculty of Science

Compact SANS Construction Proposal

Instrument Front-end

Kaspar Hewitt Klenø

Niels Bohr Institute



The Instrument Team

Kaspar Hewitt Klenø	Niels Bohr Institute
Søren Kynde	Niels Bohr Institute
Gergely Nagy	Paul Scherrer Institute
Nicholas Skar-Gislinge	Niels Bohr Institute
Kell Mortensen	Niels Bohr Institute
Kim Lefmann	Niels Bohr Institute
Joachim Kolhbrecher	Paul Scherrer Institute
Lise Arleth	Niels Bohr Institute



Concept

The instrument is designed to be a high flux broad q-range SANS intended primarily for biological samples.


This placed the following requirements on the beam delivered by the instrument front end:

- Beam divergence: up to $\pm 0.7^\circ$
- Beamspot radius: 4 mm
- No frame overlap.
- Wavelength range of 3-18.9 Å.
- Good cutoff of lower wavelengths.
- Relaxed wavelength resolution.

Sample position was fixed at 16 m from the source, as a compromise between bandwidth, wavelength resolution, and background reduction.




Overview




Guide

Width & height: 2 cm.
Start and end: 2-11.9 m.
Radius of curvature: 280 m.
Twice out of LoS.
Minimum shielding: 13.6 m.
Outer wall: $m=2.4$.
Inner wall: $m=1$.
Top and bottom: $m=2.1$.




Guide





Brilliance Transfer

Liouville transmission: 1 m coll. & $\pm 0.1^\circ$ div.




Brilliance Transfer



Collimation Section

Collimation length: 1-4 m.
1st slit radius: 8 mm.
2nd slit radius: 4 mm.
Extra guide segments when
not using the full 4 m
collimation length.

Divergence at sample




Resolution and Flux

L_1 (m)	$d\lambda/\lambda$ at 4Å	$d\lambda/\lambda$ at 8Å	$d\lambda/\lambda$ at 12Å	flux (n/s/cm ²)
1	17%	8%	6%	12×10^8
2	16%	8%	5%	3.3×10^8
4	14%	7%	5%	0.88×10^8
1 (C)	9%	8%	6%	4.7×10^8
2 (C)	8%	8%	5%	1.3×10^8
4 (C)	7%	7%	5%	0.34×10^8



Chopper System



Chopper System

Chopper Positions:

Bandwidth choppers: 6.5 & 9.5 m.

Resolution choppers: 7, 10, & 12 m.

